

## Release Note 3.45.4.2958



### 1 PURPOSE

This document describes the Release Notes for the Spirare Software, including Changes made to the SW, bug fixes, release version and date and known issues.

### 2 DOCUMENT CONTENT

#### 2.1 Version

The version number for this release is 3.45.4.2958.

New version is release due to meet the new ndd sensors.

#### 2.2 Changes

##### Spirometry

- Spirare was unable to recognize the serial numbers from the latest batch of Spirometry sensors received from the manufacturer.

##### ABPM

- An update to the FDTI driver for ABPM resulted in Spirare being unable to detect the COM port number.

#### 2.3 Bug fixes

N/A

#### 2.4 Known issues

N/A

#### 2.5 Release date

Release date for this version is 11-03-2024.