

## Release Note 3.46.2.2961



### 1 PURPOSE

This document describes the Release Notes for the Spirare PC Client Software. It includes changes to the software, bug fixes, release version and date, and known issues.

### 2 DOCUMENT CONTENT

#### 2.1 Version

The version number for this release is 3.46.2.2961.

#### 2.2 Changes

##### Various

- PC client does support Spirare web-journal's FHIR servers that generate UID-based resource IDs.

#### 2.3 Bug fixes

##### Spirometry

- Fixed an issue where the software fails to find the baseline when using the spirometry device model SPS320.

#### 2.4 Known issues

- The Child incentive feature is unstable and may sometimes require a program restart.

#### 2.5 Release date

Release date for this version is 26-04-2024.